

Total protein FS

Application for serum and plasma samples

This application was set up and evaluated by DiaSys. It is based on the standard equipment at that time and does not apply to any equipment modifications undertaken by unqualified personnel

| Identification | |
|-------------------------------------|-----|
| This method is usable for analysis: | Yes |
| Twin reaction: | No |
| Name: | TP |
| Shortcut: | |
| Reagent barcode reference: | 050 |
| Host reference: | |

| Technic | |
|----------------------------------|------------|
| Type: | End point |
| First reagent:[µL] | 180 |
| Blank reagent | Yes |
| Sensitive to light | |
| Second reagent:[µL] | 45 |
| Blank reagent | Yes |
| Sensitive to light | |
| Main wavelength:[nm] | 546 |
| Secondary wavelength:[nm] | - |
| Polychromatic factor: | - |
| 1 st reading time [min:sec] | (04:24) |
| Last reading time [min:sec] | 10:00 |
| Reaction way: | Increasing |
| Linear Kinetics | |
| Linearity: Maximum deviation [%] | |
| Fixed Time Kinetics | |
| Endpoint | |
| Prozone Limit [%] | |

| Reagents | |
|----------|--|
| Decimals | |
| Units | |

| Sample | |
|--------------------------------------|------------------|
| Diluent | DIL A (NaCl) |
| Hemolysis: | |
| Agent [µL] | 0 (no hemolysis) |
| Cleaner | |
| Sample [µL] | 0 |
| Technical limits | |
| Concentration technical limits-Lower | 0.06 |
| Concentration technical limits-Upper | 14.0 |
| SERUM | |
| Normal volume [µL] | 5 |
| Normal dilution (factor) | 1 |
| Below normal volume [µL] | 9 |
| Below normal dilution (factor) | 1 |
| Above normal volume [µL] | 5 |
| Above normal dilution (factor) | 6 |
| URIN | |
| Normal volume [µL] | 5 |
| Normal dilution (factor) | 1 |
| Below normal volume [µL] | 9 |
| Below normal dilution (factor) | 1 |
| Above normal volume [µL] | 5 |
| Above normal dilution (factor) | 6 |
| PLASMA | |
| Normal volume [µL] | 5 |
| Normal dilution (factor) | 1 |
| Below normal volume [µL] | 9 |
| Below normal dilution (factor) | 1 |
| Above normal volume [µL] | 5 |
| Above normal dilution (factor) | 6 |
| CSF | |
| Normal volume [µL] | 5 |
| Normal dilution (factor) | 1 |
| Below normal volume [µL] | 9 |
| Below normal dilution (factor) | 1 |
| Above normal volume [µL] | 5 |
| Above normal dilution (factor) | 6 |
| Whole blood | |
| Normal volume [µL] | 5 |
| Normal dilution (factor) | 1 |
| Below normal volume [µL] | 9 |
| Below normal dilution (factor) | 1 |
| Above normal volume [µL] | 5 |
| Above normal dilution (factor) | 6 |

| Results | |
|---------------------------|-------|
| Decimals | 2 |
| Units | g/dL |
| Correlation factor-Offset | 0.000 |
| Correlation factor-Slope | 1.000 |

| Range | |
|-------------|-------------|
| Gender | All |
| Age | |
| SERUM | >=6.6 <=8.8 |
| URINE | |
| PLASMA | >=6.6 <=8.8 |
| CSF | |
| Whole blood | |
| Gender | |
| Age | |
| SERUM | |
| URINE | |
| PLASMA | |
| CSF | |
| Whole blood | |

| Contaminants | |
|---------------|------------------------------|
| Contaminant 1 | Please refer to r910 Carryov |
| Wash with | |
| Cycle | |
| Volume [µL] | |
| Contaminant 2 | |
| Wash with | |
| Cycle | |
| Volume [µL] | |
| Contaminant 3 | |
| Wash with | |
| Cycle | |
| Volume [µL] | |

| Calibrators details | |
|---------------------|-----------------------|
| Calibrator list | Concentration |
| Cal. 1/Blank | 0 |
| Cal. 2 | * |
| Cal. 3 | |
| Cal. 4 | |
| Cal. 5 | |
| Cal. 6 | |
| | Max delta abs. |
| Cal. 1 | 0.002 |
| Cal. 2 | 0.020 |
| Cal. 3 | |
| Cal. 4 | |
| Cal. 5 | |
| Cal. 6 | |
| Drift limit [%] | 0.8 |
| Calculations | |
| Model | X |
| Degree | 1 |

* Enter calibrator value