PROLINE Rheumatoid factor FS

Informasi Kemasan

No. Katalog	Isi per Kit
1 7022 99 88 921	R1 4 x 18 mL + R2 4 x 4 mL
1 7022 99 88 021	R1 5 x 25 mL + R2 1 x 25 mL

Tujuan Penggunaan

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Pustaka

1. Winchester RJ. Characterization of IgG complexes in patients with rheumatoid arthritis. *Ann N Y Acad Sci* 1975; 256: 73-81.
2. Moore TL, Dorner RW. Rheumatoid factors. *Clin Biochem* 1993; 26: 75-84
3. Shmerling RH, Delbanco TL. The rheumatoid factor: an analysis of clinical utility. *Am J Med* 1991; 91: 528-34
4. Mannik M, Rheumatoid factors in the pathogenesis of rheumatoid arthritis. *J Rheumatol Suppl* 1992; 32: 46-9.
5. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. *ClinChemLabMed* 2007;45(9):1240-1243.
6. Gunder WG, Zawta B et al. *The quality of Diagnostic Samples. 1st ed. Darmstadt: Git Verlag, 2001:42-3*
7. Young DS. *Effects of Drugs on clinical Laboratory Tests - Drugs Disease, Herbs & Natural Products*, <https://clinf.wiley.com/aaccweb/aacc/>. accessed on February 2021. Published by AACC Press and John Wiley and Sons, Inc.
8. Ulvestad E, Kanestrom A, Madland TM, Thomassen E, Haga HJ. *Clinical utility of diagnostic tests for rheumatoid factor. Scandinavian Journal of Rheumatology* 2001; 30: 87-91
9. Ulvestad E, Kanestrom A, Madland TM, Thomassen E, Haga HJ. *Clinical utility of diagnostic tests for rheumatoid factor. Scandinavian Journal of Rheumatology* 2001; 30: 87-91

DiaSys Diagnostic Systems GmbH
Alte Strasse 9 65558 Holzheim
Germany

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PT Prodia Distribusi Medika
Grha Prodia Utama Lantai 3 Unit 3C
Jl.Raden Saleh No.53, Jakarta Pusat, Indonesia 10330