# PROLINE<sup>b</sup> Cystatin C FS

## Informasi Kemasan

No. Katalog	Isi per Kit
1 7158 99 10 966	R1 2 x 13 mL + R2 2 x 5 mL
1 7158 99 10 921	R1 4 x 20 mL + R2 4 x 8 mL

## Tujuan Penggunaan

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\*\* Menurut dokumen NCCLS EP17-A, vol. 24, no. 34

- Young DS. *Effects of Drugs on Clinical Laboratory Tests*. 5<sup>th</sup> ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
- Young DS. *Effects on Clinical Laboratory Tests – Drugs Disease, Herbs & Natural Products*, <https://clinf.wiley.com/aaccweb/aacc/>, accessed in August 2021. Published by AACC Press and John Wiley and Sons, Inc.
- Bökenkamp A, Grabensee A, Stoffel-Wagner B, Hasan C, Henne T, Offner G, Lentze MJ. The  $\beta$ 2-microglobulin/cystatin C ratio: a potential marker of post-transplant lymphoproliferative disease. *Clin Nephrol* 2001; 58: 417-22.
- Randers E, Krue S, Erlandsen EJ, Danielsen H, Hansen H, Hansen LG. Reference interval for serum cystatin C in children. *ClinChem* 1999; 45:1856-60.
- Finney H, Newman DJ, Price CP. Adult reference ranges for serum cystatin C, creatinine and predicted creatinine clearance. *Ann Clin Biochem* 2000; 37: 49-59.

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## Pustaka

- Erland J.E., Randers E., Kristensen J.H. Reference intervals for serum cystatin C and serum creatinine in adults. *Clin Chem Lab Med* 1998; 36(6):393-397.
- Lamb E., Newman DJ, Price CP. *Kidney function tests*. In: Burtis CA, Ashwood ER, Bruns DE, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*. 4th edition St. Louis Missouri: Elsevier Saunders; 2006; p. 823-835.
- Kyhse-Andersen, Schmidt C., Nordin G. et al. Serum cystatin C, determined by a rapid, automated particle-enhanced turbidimetric method, is a better marker than serum creatinine for glomerular filtration rate. *ClinChem* 1994; 40(10):1921-6.
- Le Bricon T., Leblanc I et al. Evaluation of renal function in intensive care: plasma cystatin C vs. creatinine and derived glomerular filtration rates *Clin Chem Lab Med* 2005; 43(9):953-957.
- Le Bricon T., Thervet E., Benlakehal M. et al. Changes in Plasma cystatin C after renal transplantation and acute rejection in adults. *Clin Chem* 1999; 45(12):2243-9
- Ustundag Y., Samsar U. et al. Analysis of glomerular filtration rate, serum cystatin C levels, and resistive index values in cirrhosis patients. *Clin Chem Lab Med* 2007; 45(7):890-94.
- Stabuc B., Vrhovec L. et al. Improved prediction of decreased creatinine clearance by serum cystatin C: Use in cancer patients before and during chemotherapy. *Clin Chem* 2000; 46(2):193-7.
- Pucci L., Triscornia S., Lucchesi D. et al. Cystatin C and estimates of renal function: searching for a better measure of kidney function in diabetic patients. *Clin Chem* 2007; 53(3):480-8.
- Strevens H, Wide-Svensson D, Grubb A. Serum cystatin C is a better marker for preeclampsia than serum creatinine or serum urate. *Scand J Clin Lab Invest* 2001; 61:575-80.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. *ClinChemLabMed* 2007;45(9):1240-1243.
- Guder WG, Zawta B. et al. *The Quality of Diagnostic Samples*. 1st ed. Darmstadt: GIT Verlag; 2001; p. 24-5.
- Wiesli P., Schwegler B. et al. Serum cystatin C is sensitive to small changes in thyroid function. *Clinica Chimica Acta* 2003; 338: 87-90.

DiaSys Diagnostic Systems GmbH  
Germany

### Diimpor dan didistribusikan oleh:

PT Prodia Diagnostic Line  
Cikarang, Indonesia