# PROLINE<sup>b</sup> Procalcitonin FS

## Informasi Kemasan

No. Katalog                      Isi per Kit  
1 7318 99 10 930                R1: 2 x 18 mL + R2: 2 x 6 mL

## Tujuan Penggunaan

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## Pustaka

- Rudd KE et al. Global, regional, and national sepsis incidence and mortality, 1990–2017: analysis for the Global Burden of Disease Study. *The Lancet* 2020; 395 (10219):200-211.
- Singer M et al. The Third International Consensus Definitions for Sepsis and Septic Shock (Sepsis-3). *JAMA* 2016; 315(8): 801-810.
- Fleischmann C et al. Assessment of global incidence and mortality of hospital-treated sepsis. *Current estimates and limitations*. *Am J Respir Crit Care Med*. 2016; 193(3):259-272.
- Maruna P, Nedeljković K and Gürlich R. Physiology and genetics of procalcitonin. *Physiol Res*. 2000; 49(Suppl 1): S57–S61.
- Christ-Crain M, Müller B. Procalcitonin in bacterial infections-hype, hope, more or less? *Swiss Med Weekly*. 2005; 135:451-460.
- Becker KI et al. Procalcitonin in sepsis and systemic inflammation: a harmful biomarker and a therapeutic target. *British journal of pharmacology* 2010; 159(2): 253-264.
- Becker KI et al. Procalcitonin and the calcitonin gene family of peptides in inflammation, infection, and sepsis: a journey from calcitonin back to its precursors. *The Journal of Clinical Endocrinology & Metabolism*. 2004; 89(4): 1512-1525.
- Müller B et al. Ubiquitous expression of the calcitonin-i gene in multiple tissues in response to sepsis. *J Clin Endocrinol Metab* 2001; 86(1): 396-404.
- Bartoletti, Michele, et al. Procalcitonin-guided antibiotic therapy: an expert consensus. *Clin Chem Lab Med*. 2018;56:1223-1229.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. *ClinChemLabMed* 2007; 45(9): 1240-1243.
- Guzdys V et al. Method Verification Shows a Negative Bias between 2 Procalcitonin Methods at Medical Decision Concentrations. *The journal of applied laboratory medicine* 2019; 4(1): 69-77.
- Meisner M. Procalcitonin-influence of temperature, storage, anticoagulation and arterial or venous asseration of blood samples on procalcitonin concentrations. *Clinical Chemistry and Laboratory Medicine* 1997; 35(8): 597-602.
- Young DS. *Effects of Drugs on Clinical Laboratory Tests*. 5<sup>th</sup> ed. Volume 1 and 52. Washington, DC: The American Association for Clinical Chemistry Press 2000.
- Young DS. *Effects on Clinical Laboratory Tests - Drugs, Disease, Herbs & Natural Products*, <https://clinfx.wiley.com/aaccweb/aacc/>, accessed in June 2021. Published by AACC Press and John Wiley and Sons, Inc
- Harbarth S et al. Diagnostic value of procalcitonin, interleukin-6, and interleukin-8 in critically ill patients admitted with suspected sepsis. *Am J Respir Crit Care Med* 2001; 164:396-402.
- Meisner M. Procalcitonin - *Biochemie und klinische Diagnostik*. 1. Auflage Bremen: UNI-MED- Verlag 2010.
- Chiesa C et al. Reliability of procalcitonin concentrations for the diagnosis of sepsis in critically ill neonates. *Clinical infectious diseases* 1998; 26(3): 664-672.
- Chiesa C et al. C-reactive protein, interleukin-6, and procalcitonin in the immediate postnatal period: influence of illness severity, risk status, antenatal and perinatal complications, and infection. *Clinical chemistry* 2003; 49(1):60-68

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DiaSys Diagnostic Systems GmbH