Proline CRP U-hs

Informasi Kemasan

No. Katalog	Isi per Kit
1 7045 99 88 920	R1 4 x 11 mL + R2 4 x 11 mL
1 7045 99 88 930	R1 3 x 20 mL + R2 3 x 20 mL

Tujuan Penggunaan

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Pustaka

1. Thompson D, Milford-Ward A, Whicher JT. The value of acute phase protein measurements in clinical practice. *Ann Clin Biochem* 1992;29:123-31.
2. Gabay C, Kushner I. Acute-phase proteins and other systemic responses to inflammation. *N Engl J Med* 1999;340:448-54.
3. Hansson LO, Lindquist L. C-reactive protein: its role in the diagnosis and follow-up of infectious diseases. *Curr Opin Infect Diseases* 1997; 10:196–201.
4. Sipe JD. Acute-phase proteins in osteoarthritis. *Semin Arthritis Rheum* 1995;25: 75–86.
5. Guder WG, Zawta B et al. *The Quality of Diagnostic Samples*. 1st ed. Darmstadt: GIT Verlag; 2001. p. 24-5.
6. Dati F, Schumann G, Thomas L, Aguzzi F, Baudner S, Bienvenu J et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP reference material (CRM 470). *Eur J Clin Chem Clin Biochem* 1996;34:517-20.
7. Schlebusch H, Liappis N, Klein G. High sensitive CRP and creatinine: reference intervals from infancy to childhood Poster presented at AACC/CSCC; July/August 2001, Chicago, Illinois.
8. Young DS. *Effects of Drugs on Clinical Laboratory Tests*. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
9. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. *ClinChemLabMed* 2007;45(9):1240–1243.

DiaSys Diagnostic Systems GmbH

Alte Strasse 9 65558 Holzheim

Germany

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PT Prodia Distribusi Medika

Grha Prodia Utama Lantai 3 Unit 3C

Jl. Raden Saleh No. 53, Jakarta 10330 - Indonesia