Proline Pancreatic Amylase CC FS

Informasi Kemasan

No. Katalog	Isi per Kit
1 0551 99 88 021	R1 5 x 20 mL + R2 1 x 25 mL
1 0551 99 88 921	R1 4 x 20 mL + R2 4 x 5 mL

Tujuan Penggunaan

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Pustaka

1. Lorentz K. α -Amylase. In: Thomas L, editor. *Clinical laboratory diagnostics*. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 192-202.
2. Moss DW, Henderson AR. *Digestive enzymes of pancreatic origin*. In: Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 689-98.
3. Gerber M, Naujoks K, Lenz H, Wulff K. A monoclonal antibody that specifically inhibits human salivary α -amylase. *Clin Chem* 1987; 33: 1158-62.
4. Kruse-Jarres JD, Kaiser C, Hafkenschied JC, Hohenwallner W, Stein W., Bohner J et al. Evaluation of a new α -amylase assay using 4,6-ethylidene-(G7)-1-4-nitrophenyl-(G1)- α , D-maltoheptaoside as substrate. *J Clin Chem Biochem* 1989; 27: 103-13.
5. Tietz NW, Burlina A, Gerhardt W, Junge W, Mafferttheimer P, Mural T et al. Multicenter evaluation of a specific pancreatic isoamylase assay based on a double monoclonal-antibody technique. *Clin Chem* 1988; 34: 2096-102.
6. Guder WG, Zawta B et al. *The Quality of Diagnostic Samples*. 1st ed. Darmstadt: GIT Verlag; 2001. p. 16-17.
7. Young DS. *Effects of Drugs on Clinical Laboratory Tests*. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press, 2000.
8. Junge W, Wortmann W, Wilke B, Waldenstroem J et al. Development and evaluation of assays for determination of total and pancreatic amylase at 37 °C according to the principle recommended by the IFCC. *Clin Biochem* 2001; 34: 607-15.
9. Junge W, Troge B, Klein G, Poppe W, Gerber M. Evaluation of a new assay for pancreatic amylase: Performance characteristics and estimation of reference interval. *Clin Biochem* 1989; 22: 109-14.
10. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: Mechanisms, detection and prevention. *Clin Chem Lab Med* 2007; 45(9): 1240-1243.

  **DiaSys Diagnostic Systems GmbH**
Alte Strasse 9 65558 Holzheim
Germany

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PT Prodia Distribusi Medika

Grha Prodia Utama Lantai 3 Unit 3C

Jl. Raden Saleh No. 53, Jakarta 10330 - Indonesia