Proline Pancreatic Amylase CC FS

Informasi Kemasan

No. Katalog	Isi per Kit
1 0551 99 88 021	R1 5 x 20 mL + R2 1 x 25 mL
1 0551 99 88 921	R1 4 x 20 mL + R2 4 x 5 mL

Tujuan Penggunaan

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Pustaka

- Moss DW, Henderson AR. Digestive enzymes of pancreatic origin. In: Burtis CA, Ashwood ER, editors. *Tietz Textbook of Clinical Chemistry*. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 689-98.
- Thomas L. *Clinical Laboratory Diagnostics* [Internet]; 2020 [cited 2023 Jun 19]. Available from: <https://www.clinical-laboratory-diagnostics-2020.com/>
- Gerber M, Naujoks K, Lenz H, Wulff K. A monoclonal antibody that specifically inhibits human salivary alpha-amylase. *Clin Chem* 1987; 33: 1158-62.
- Kruse-Jarres JD, Kaiser C, Hafkenscheid JC, Hohenwallner W, Stein W., Bohner J et al. Evaluation of a new alpha-amylase assay using 4,6- ethylidene-(G7)-1-4-nitrophenyl-(G1)-alpha, D-maltoheptaoside as substrate. *J Clin Chem Biochem* 1989; 27: 103-13.
- Kruse-Jarres JD, Burlina A, Gerhardt W, Junge W, Maffetherheimer P, Mural T et al. Multicenter evaluation of a specific pancreatic isoamylase assay based on a double monoclonal-antibody technique. *Clin Chem* 1988; 34: 2096–102
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: Mechanisms, detection and prevention. *Clin Chem Lab Med* 2007; 45(9): 1240–1243
- Guder WG, Zawta Bet al. *The Quality of Diagnostic Samples*. 1st ed. Darmstadt: GIT Verlag; 2001. p. 16-17. 50-51.
- Young DS. *Effects of Drugs on Clinical Laboratory Tests*. 5th. ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press, 2000.
- Young DS. *Effects on Clinical Laboratory Tests - Drugs Disease, Herbs & Natural Products*, <https://clinfx.wiley.com/aaccweb/aacc/>, accessed in February 2024. Published by AACC Press and John Wiley and Sons, Inc.
- Junge W, Wortmann W, Wilke B, Waldenstroem J et al. Development and evaluation of assays for determination of total and pancreatic amylase at 37 °C according to the principle recommended by the IFCC. *Clin Biochem* 2001; 34: 607-15.

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DiaSys Diagnostic Systems GmbH