

Updated sheets: [www.proline.co.id/document/](http://www.proline.co.id/document/)

Value sheet version: 2

| Short Name / Analyte   | Method               | Reagent ID | Traceability   | Value | Unit   |
|--|----------------------|------------|--|-------|--------|
| <b>HDL-C</b><br>HDL cholesterol  | Immunoinhibition     | 1 3521     | NIST SRM 1951 Level 2  | 64.5  | mg/dL  |
|  |                      |            |  | 1.67  | mmol/L |
|  |                      |            |  | 0.645 | g/L    |
|  | Direct, colorimetric | 1 3561     | Commercially available assay traceable to CDC reference method | 60.9  | mg/dL  |
|  |                      |            |  | 1.57  | mmol/L |
|  |                      |            |  | 0.609 | g/L    |
| <b>LDL-C</b><br>LDL cholesterol  | Select               | 1 4121     | NIST SRM 1951b Level 2   | 140   | mg/dL  |
|  |                      |            |  | 3.61  | mmol/L |
|  |                      |            |  | 1.40  | g/L    |
|  | Direct, colorimetric | 1 4131     | NIST SRM 1951c Level 2   | 149   | mg/dL  |
|  |                      |            |  | 3.85  | mmol/L |
|  |                      |            |  | 1.49  | g/L    |
| <b>Lp-PLA<sub>2</sub></b><br>Lipoprotein associated phospholipase A <sub>2</sub> | Colorimetric         | 1 7181     | ε 4-nitrophenol  | 426   | U/L    |
|  |                      |            |  | 7.11  | μkat/L |
| <b>NEFA</b><br>Non esterified fatty acids  | Enzymatic            | 1 5781     | Internal Reference   | 51.6  | mg/dL  |
|  |                      |            |  | 1.83  | mmol/L |
|  |                      |            |  | 0.516 | g/L    |
| <b>PL</b><br>Phospholipids   | Enzymatic            | 1 5741     | Internal Standard  | 318   | mg/dL  |
|  |                      |            |  | 4.10  | mmol/L |
|  |                      |            |  | 3.18  | g/L    |

\*This value sheet is short version. For other methods, please ask your distributors or email to [support@proline.co.id](mailto:support@proline.co.id)

**Relabeled and distributed by:**

 PT Prodia Diagnostic Line  
 Kawasan Industri Jababeka III  
 Jl. Tekno 1 Blok C 2 D-E-F  
 Cikarang, Jawa Barat 17530 - Indonesia