

Short Name / Analyte	Method*	Reagent Cat. No.	Traceability	Value	Unit
ALB Albumin	Bromcresol green	1 0220	ERM-DA 470	3.76 37.6	g/dL g/L
AP Alkaline Phosphatase	IFCC 37°C	1 0441	IFCC recommended formulation 37°C	234 3.91	U/L μkat/L
ALT/GPT Alanine aminotransferase	IFCC	1 2701 +	Original IFCC formulation	96.0	U/L
	with pyridoxal phosphate 37°C	2 5010	37°C	1.60	μkat/L
	w/o pyridoxal phosphate 37°C	1 2701	Modified IFCC formulation ε NADH, 37°C	96.8 1.61	U/L μkat/L
AST/GOT Aspartate aminotransferase	IFCC	1 2601 +	Original IFCC formulation	96.7	U/L
	with pyridoxal phosphate 37°C	2 5010	37°C	1.61	μkat/L
	w/o pyridoxal phosphate 37°C	1 2601	Modified IFCC formulation ε NADH, 37°C	97.8 1.63	U/L μkat/L
AMY Amylase	EPS G-7	1 0501	Original IFCC formulation 37°C	174 2.90	U/L μkat/L
			Modified IFCC formulation ε p-Nitrophenol, 37°C	157 2.62	U/L μkat/L
DBIL Direct Bilirubin	DCA	1 0821	Jendrassik-Grof	2.38 40.7 23.8	mg/dL μmol/L mg/L
TBIL Total Bilirubin	DCA	1 0811	SRM 916a	3.95 67.5 39.5	mg/dL μmol/L mg/L
BUN Urea N	Urease UV	1 3101	SRM 909b Level 1	49.1 17.5 0.491	mg/dL mmol/L g/L
CAAS Calcium	Arsenazo	1 1130	AAS	2.60 5.20 10.4	mmol/L mval/L mg/dL
CA Calcium	Phosphonazo	1 1181	AAS	2.80 5.60 11.2	mmol/L mval/L mg/dL
CHOL Cholesterol	CHOD-PAP	1 1300	GC-IDMS	153 3.96 1.53	mg/dL mmol/L g/L
CHE Cholinesterase	Butyrylthiocholin	1 1401	DGKC recommended formulation 37°C	5137 85.6	U/L μkat/L
CK Creatine kinase	DGKC opt.	1 1601	Original IFCC formulation 37°C	370	U/L
	IFCC 37°C			6.16	μkat/L
CREAP Creatinine	Enzymatic PAP	1 1759	GC-IDMS	3.87 342 38.7	mg/dL μmol/L mg/L
CREA Creatinine	Jaffé, w/o compensation	1 1711	GC-IDMS	3.87 342 38.7	mg/dL μmol/L mg/L
	Jaffé, with compensation	1 1711	GC-IDMS	4.32 382 43.2	mg/dL μmol/L mg/L
FE Iron	Ferene	1 1911	SRM 682	192 34.3 1.92	μg/dL μmol/L mg/L
	Ferene respons@920	1 1911	SRM 682	198 35.4	μg/dL μmol/L
				1.98	mg/L

Short Name / Analyte	Method*	Reagent Cat. No.	Traceability	Value	Unit
GGT γ-Glutamyltransferase	Szasz liquid stand. Szasz	1 2801	Original formulation	83.9	U/L
			Persijn/v.d.Slik 1976, 37°C	1.40	μkat/L
	Szasz liquid stand. IFCC	1 2801	Original IFCC formulation	94.9	U/L
			37°C	1.58	μkat/L
GLDH Glutamate dehydrogenase	DGKC 1970 37°C	1 2411	DGKC recommended formulation	27.4	U/L
				0.456	μkat/L
GLUC/GLUHK Glucose	GOD-PAP (10 min.) HK/G6P-DH	1 2500	ID-MS	194	mg/dL
		1 2511		10.7	mmol/L
				1.94	g/L
LACT Lactate	LDH UV	1 4001	Primary reference material	27.8	mg/dL
				3.08	mmol/L
				278	mg/L
LDH21 Lactate dehydrogenase	IFCC, modified	1 4251	Original IFCC formulation	262	U/L
				4.37	μkat/L
LPS Lipase	Colorimetric	1 4321	Molar absorbance coefficient	102	U/L
				1.69	μkat/L
MG Magnesium	Xylidyl blue	1 4610	AAS	1.15	mmol/L
				2.30	mval/L
				2.79	mg/dL
PO4 Phosphorus	Molybdate UV	1 5211	Internal Standard	5.81	mg/dL
				1.87	mmol/L
				58.1	mg/L
TP Total Protein	Biuret, with sample blank	1 2311	Modified Biuret method	5.53	g/dL
				55.3	g/L
	Biuret, w/o sample blank	1 2311	Modified Biuret method	6.06	g/dL
				60.6	g/L
TRIG Triglycerides	GPO-PAP	1 5710	GC-IDMS	135	mg/dL
				1.52	mmol/L
				1.35	g/L
UIBC Unsat. Iron Binding Cap.	Ferene / Ferrozine	1 1921	Internal Standard	216	μg/dL
				38.8	μmol/L
				2.16	mg/L
UREA Urea	Urease UV	1 3101	SRM 909b Level 1	105	mg/dL
				17.5	mmol/L
				1.05	g/L
UA Uric Acid	Uricase - TOOS	1 3001	GC ID-MS	5.15	mg/dL
				306	μmol/L
				51.5	mg/L
UATBH Uric Acid	Uricase - TBHBA	1 3021	GC ID-MS	5.15	mg/dL
				306	μmol/L
				51.5	mg/L

*This value sheet is short version. For other methods, please ask your distributors or email to technical.support@proline.co.id

IVD  **DiaSys Diagnostic Systems GmbH**
Alte Strasse 9 65558 Holzheim
Germany

Repacked by
PT Prodia Diagnostic Line
Kawasan Industri Jababeka III
Jl. Tekno 1 Blok C 2 D-E-F
Cikarang, Kab.Bekasi, Jawa Barat 17530 - Indonesia

Distributed by
PT Prodia Distribusi Medika
Grha Prodia Utama Lantai 3 Unit 3C
Jl. Raden Saleh No. 53, Jakarta 10330
Indonesia